
ENG/RUS

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DIRECTIONS

- ☐ All the time you are in the laboratory you must wear safety spectacles or optical spectacles, in case the latter are needed. Pipettes have to be filled with a pipette bulb. Eating in the laboratory is strictly prohibited.
- ☐ Participants have to follow the safety rules, to be kind to each other and to keep their labware and working place in order. Do not be shy to ask the laboratory assistant about the safety rules.
- ☐ You may begin working only after the “START” signal. You have 5 hours to complete the practical examination and to fill in the answer sheets. 15 minutes before the end of this period you will be warned. Working has to be stopped right after the “STOP” signal. Being 5 minutes late will result in receiving 0 points for the practical examination.
- ☐ Write your name and code (indicated in the working place) in the corresponding boxes on the answer sheets.
- ☐ All results have to be written in the corresponding boxes on the answer sheets. Results written elsewhere will not be graded. If you need some paper for calculations or an additional answer sheet, ask the laboratory assistant.
- ☐ When you have finished working, place all answer sheets into the provided envelope. Only the sheets in the envelope will be graded.
- ☐ Do not leave the room prior to receiving permission.
- ☐ You have to use only the provided labware.
- ☐ The number of significant figures in the numerical answers has to correspond to the laws of experimental errors and rounding up. Failing to do the calculations correctly will result in penalties, even if your technique is impeccable.

Collection of used and spilt chemicals and broken glassware

- ☐ Organic filtrates, washings and any other waste have to be put into the corresponding vessels or bottles.
- ☐ For disposing of waste and chemical waste use the corresponding containers.
- ☐ Broken glassware and split glass has to be disposed of in the waste basket. For breaking glassware or asking an additional sample you receive a 2 point penalty.

Periodic Table of Elements

Atomic weights based on $^{12}\text{C} = 12$
(Numbers) = most stable isotope

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
Labware:

technical balance	1
50 mL round bottom flask	1
heating plate	1
reflux condenser	1
boiling stones	2-3
Bunsen flask	1
Buchner funnel	1
Petri dish	1
spatula	1
50 mL cylinder	1
10 mL measuring pipette	1
laboratory stand	1

In the free area to the right draw schematically the reaction assembly.

Reagents:

Ethanol

Reaction assembly

Safety information:

R11	Highly flammable
S7	Keep container tightly closed
S16	Keep away from sources of ignition - No smoking

Ethyl acetoacetate

Safety information:

R36/37/38	Irritating to eyes, respiratory system and skin
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S37/39	Wear suitable gloves and eye/face protection

Ammonium acetate

Safety information:

S24/25	Avoid contact with skin and eyes
--------	----------------------------------

Urotropine

Safety information:

R11	Highly flammable
R42/43	May cause sensitization by inhalation and skin contact
S16	Keep away from sources of ignition - No smoking
S22	Do not breathe dust
S24	Avoid contact with skin
S37	Wear suitable gloves

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Reagents and safety information:

hydrogen peroxide, H₂O₂; approximately 3% solution

Safety information:

R22	Harmful if swallowed
R37/38	Irritating to respiratory system and skin
R41	Risks of serious damage to eyes
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28A	After contact with skin, wash immediately with plenty of water
S36/37/38	Wear suitable protective clothing, gloves and eye/face protection
S45	In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible)

potassium permanganate, KMnO₄; 0.020 M solution

Safety information:

R22	Harmful if swallowed
R50/53	Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment
R8	Contact with combustible material may cause fire
S60	This material and its container must be disposed of as hazardous waste
S61	Avoid release to the environment. Refer to special instructions / Safety data sheets

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sulfuric acid, H_2SO_4 ; 1 M solution + Mn^{2+} ions

Safety information:

R35	Causes severe burns
R51/53	Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S30	Never add water to this product
S45	In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible)
S61	Avoid release to the environment. Refer to special instructions/safety data sheet

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